


TA-01 Mouse Anti-Myc mAb

Immunogen	AEEQKLISEEDLLRKRREQLKHKLE conjugated to KLH, corresponding to C terminal amino acids 408-432 of human c-Myc
Product name	Mouse anti-Myc monoclonal antibody
Clone	9E10
Isotype	IgG1
Amount	0.1 ml
Reactivity	Myc Tag
Storage/Stability	Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Avoid repeated freeze-thaws.
Application	WB
Suggested dilutions	1:1000-1:2000

Experimental data



Immunoblot analysis of four different Myc-DDK tagged over-expression lysates and one Empty vector transfected 293T cells with Anti-Myc at 1:1000.

Notes: This product is to be used for laboratory only. Not for diagnostic use.

