




TA-02 Mouse Anti-His mAb

Immunogen	Synthetic peptide (HHHHHH) coupled to KLH
Product name	Mouse anti-His monoclonal antibody
Clone	OTI2B5
Isotype	IgG1
Amount	0.1 ml
Reactivity	His tag
Storage/Stability	Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Avoid repeated freeze-thaws.
Application	WB
Suggested dilutions	1:1000-1:2000

Experimental data



Anti-His Immunoblot analysis of six different HEK293T over-expression lysates containing different His tagged recombinant proteins, and two his tagged bacterial lysates produced in E. coli at 1:1000 dilution.

Notes: This product is to be used for laboratory only. Not for diagnostic use.



北京中杉金桥生物技术有限公司 Http://www.zsbio.com 地址:北京市西城区广外小马厂东里 10 号智恒泰大厦 2 层 电话: 010-63470285 400-810-9781 E-mail: sales@zsbio.com

