




## TA-05 Mouse Anti-DDK(FLAG) mAb

Immunogen	Synthetic peptide (DYKDDDDK) coupled to KLH
Product name	Mouse anti-DDK (FLAG) monoclonal antibody
Clone	OTI11C3
Isotype	IgG1
Amount	0.1 ml
Reactivity	DDK(FLAG Tag): DYKDDDDK
Storage/Stability	Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Avoid repeated freeze-thaws.
Application	WB
Suggested dilutions	1:1000-1:2000

## **Experimental data**



HEK293T cells were transfected with the pCMV6-ENTRY YPEL3 (RC208686; lane1), CYSLTR1 (RC207763; lane2), SEMA3D (RC216032; lane3), Dmp1 (RC222208; lane4) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DDK (1:1000).

\* Peptide sequence of the DDK-tag (Flag®): N-DYKDDDDK-C Flag® is a registered trademark of Sigma-Aldrich.

Notes: This product is to be used for laboratory only. Not for diagnostic use.



北京中杉金桥生物技术有限公司 Http://www.zsbio.com 地址:北京市西城区广外小马厂东里10号智恒泰大厦2层 电话: 010-63470285 400-810-9781 E-mail: sales@zsbio.com

