

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA309729

LDL Receptor (LDLR) Mouse Monoclonal Antibody [Clone ID: C7]

Product data:


Product Type:	Primary Antibodies
Clone Name:	C7
Applications:	IF, IHC
Recommend Dilution:	WB 1:1000; IF/ICC 1:100; IHC 1:200; IP 1:10-1:500
Reactivity:	Bovine, Human
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	Partially purified bovine LDL receptor [Swiss-Prot# P01131]
Formulation:	Tris-glycine, 150mM NaCl pH 7.5 and 0.05% Sodium Azide
Purification:	Protein G purified
Gene Name:	low density lipoprotein receptor
Database Link:	NP 000518 Entrez Gene 3949 Human
Synonyms:	FH; FHC; LDLCQ2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Endocytosis




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

CRIGENE LDL Receptor (LDLR) Mouse Monoclonal Antibody [Clone ID: C7] – TA309729

Product images:



Immunohistochemistry: LDL Receptor Antibody (C7) - Immunohistochemical analysis of LDL Receptor on human liver using anti-LDL Receptor.



Immunocytochemistry/Immunofluorescence: LDL Receptor Antibody (C7) - Antibody was tested in HepG2 cells with FITC (green). Nuclei were counterstained with Dapi (blue).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US