

## **Product datasheet for TA309827**

## RAD51D Mouse Monoclonal Antibody [Clone ID: 5B3/6]

## **Product data:**

**Product Type:** Primary Antibodies

Clone Name: 5B3/6
Applications: WB

Recommend Dilution: WB, IF, ICC

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** His-tagged hRAD51D, overexpressed in E. coli.

**Formulation:** PBS at physiological pH and Sodium Azide

**Purification:** Ascites

Gene Name: RAD51 paralog D

Database Link: NP 002869 Entrez Gene 5892 Human

Synonyms: BROVCA4; R51H3; RAD51L3; TRAD

**Protein Families:** Druggable Genome

**Protein Pathways:** Homologous recombination

OriGene Technologies, Inc.


Rockville, MD 20850, US Phone: +1-888-267-4436

techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

9620 Medical Center Drive, Ste 200



## **Product images:**



Rad51D detected in HEK293 lysates using a 1:1,000 dilution of anti-Rad51D (purified). Photo courtesy of B.T. Bennett & K. Knight of University of Massachusetts Medical School