

## Product datasheet for **TA309899**

### TRF1 (TERF1) Mouse Monoclonal Antibody [Clone ID: 57-6]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	57-6
Applications:	WB
Recommend Dilution:	WB, IF, ICC, ELISA
Reactivity:	Human, Mouse, Primate, Rat
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant TRF1
Formulation:	Tris-glycine, 150mM NaCl and 0.05% Sodium Azide
Concentration:	1.5mg/ml
Purification:	Protein G purified
Gene Name:	telomeric repeat binding factor 1
Database Link:	<a href="#">NP_003209 Entrez Gene 21749 Mouse</a> <a href="#">Entrez Gene 297758 Rat</a> <a href="#">Entrez Gene 7013 Human</a>
Synonyms:	hTRF1-AS; PIN2; t-TRF1; TRBF1; TRF; TRF1
Protein Families:	Transcription Factors




[View online »](#)

## Product images:



Western Blot: TRF1 Antibody (57-6) - Western blot detection of TRF1. 50 ugs of total lysate each lane.



Western Blot: TRF1 Antibody (57-6) - Detection of TRF1 in HeLa whole cell extracts