

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA309904

mtTFA (TFAM) Mouse Monoclonal Antibody [Clone ID: 18G102B2E11]

Product data:

Product Type: Primary Antibodies

Clone Name: 18G102B2E11

Applications: IF, IHC, WB

Recommend Dilution: WB, IHC, ICC

Reactivity: Human

Host: Mouse

Isotype: IgG

Clonality: Monoclonal

Immunogen: Human mtTFA [Swiss-Prot# Q00059]

Formulation: Tris-glycine, 150mM NaCl pH 7.5 and 0.05% Sodium Azide

Concentration: 1 mg/ml

Purification: Protein G purified

Gene Name: transcription factor A, mitochondrial

Database Link: NP 003192 Entrez Gene 7019 Human


Synonyms: MTTF1; MTTFA; TCF6; TCF6L1; TCF6L2; TCF6L3

Protein Families: Druggable Genome, Transcription Factors


Protein Pathways: Huntington's disease



Product images:





Western Blot: mtTFA Antibody - Western blot analysis of mtTFA expression in HeLa whole cell lysate using anti-mtTFA.



Immunohistochemistry: mtTFA Antibody - Immunohistochemical analysis of mtTFA on human testis using anti-mtTFA.







ICC: mtTFA Antibody - Immunostaining in HeLa cells. C1: DAPI staining 2 nuclei, C2: anti-mtTFA stain with Alexa Fluor 488 secondary ab (green), C3: Mitotracker Red staining mitochondria. Image courtesy of Elizabeth Wang.

Immunocytochemistry: mtTFA Antibody - mtTFA antibody was tested in HeLa cells with FITC (green). Nuclei and actin were counterstained with Dapi (blue) and Phalloidin (red).