

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA319599

IL17F Mouse Monoclonal Antibody [Clone ID: 4H1629.1]

Product data:

Product Type: Primary Antibodies

Clone Name: 4H1629.1
Applications: IHC, WB

Recommend Dilution: WB: 0.5ug/mL, IHC: 5ug/mL

Reactivity: Human Host: Mouse

Clonality: Monoclonal

Immunogen: Anti-IL-17F (MOUSE) Monoclonal Antibody was produced in mouse by repeated

immunizations with mature full length recombinant human IL-17F produced in E.coli followed

by hybridoma development.

Formulation: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Concentration: 1 mg/ml

Gene Name: interleukin 17F

Database Link: NP 443104 Entrez Gene 112744 Human

Synonyms: CANDF6; IL-17F; ML-1; ML1

Note: Anti-L-17F recognizes IL-17F (also known as Cytokine ML-1 or Interleukin-24). IL-17F is

produced and secreted by CD8+ T cells, NK cells, NKT cells and LTi cells. The main functions

of IL-17F are neutrophil recruitment and immunity to extracellular pathogen. More

importantly, IL-17F drives inflammation and auto-immunity. IL-17A and IL-17F are by far the best characterized cytokines of the IL-17 cytokine family. IL-17F dimerizes in a parallel fashion similar to nerve growth factor and other neutrophins. Its dimerization is critical to fulfill its activity. When secreted by activated T cells, IL-17F can stimulate the production of other cytokines such as IL-6, IL-8 granulocyte colony-stimulating factor and, can stimulate cartilage matrix turnover. Defects in IL17F are the cause of familial candidiasis type 6 (CANDF6). CANDF6 is a rare disorder with altered immune responses and impaired clearance of fungal infections, selective against Candida. Anti-IL-17E cytokine antibody is ideal for investigators


involved in Immunology, Signal Transduction research, Cancer and Inflammatory pathologies.

Protein Families: Druggable Genome, Secreted Protein






Product images:



Western Blot of Mouse Anti-IL-17F antibody Lane 1: human full length recombinant IL-17F protein Lane 2: mouse full length recombinant IL-17F protein Lane 3: rat full length recombinant IL-17F protein Load: 20 ng/lane Primary antibody: Anti-IL-17F antibody at 0.5ug/mL for overnight at 4°C



Immunohistochemistry of Mouse Anti-IL-17F antibody Tissue: human colon tissue Fixation: formalin-fixed, paraffin-embedded Primary antibody: isotype control (top left), Mouse Anti-IL-17F antibody (bottom left, right) at 5 ug/ml