

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA319820

# HAX1 Mouse Monoclonal Antibody [Clone ID: 9G6C6]


## **Product data:**

Product Type:	Primary Antibodies
Clone Name:	9G6C6
Applications:	WB
<b>Recommend Dilution:</b>	WB: 1 - 2 ug/mL, IF: 20 ug/mL
Reactivity:	Human, Rat
Host:	Mouse
lsotype:	lgG3
Clonality:	Monoclonal
Immunogen:	Mouse monoclonal Hax1a antibody was raised against a 15 amino acid synthetic peptide near the amino terminus of human Hax1a.
Formulation:	Hax1a Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration:	1ug/ul
Purification:	Hax1a Monoclonal Antibody is immunoaffinity chromotography purified lgG.
Gene Name:	HCLS1 associated protein X-1
Database Link:	NP_001018238 Entrez Gene 291202 RatEntrez Gene 10456 Human
Background:	Hax1a Monoclonal Antibody: The HS-1 associated protein X-1 (Hax1) was initially identified in a yeast two-hybrid assay on the basis of its ability to bind to the hemapoietic cell-specific protein 1 (HS-1). Hax1 possesses anti-apoptotic activity and is structurally related to Bcl-2 family members, including the presence of BH1- and BH2-like domains. It has recently been shown to interact with HIV viral protein R (Vpr), a protein required for viral pathogenesis of HIV and linked to T-cell apoptosis through activation of caspases 3 and 9. Other studies indicate that Hax1-mediated processing of HtrA2 (also known as Omi) by the mitochondrial protease PARL allows survival of lymphocytes and neurons when cytokines are limiting. At least four isoforms of Hax1 are known to exist. This antibody is expected to recognize the longest isoform (Hax1a) as well as the shortest.
Synonyms:	HCLSBP1; HS1BP1; SCN3



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

### **Product images:**



Western blot analysis of Hax1a in human heart tissue lysate with Hax1a antibody at (A) 1 and (B) 2 ug/mL.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US