

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA319869

MD2 (LY96) Mouse Monoclonal Antibody [Clone ID: 1A2E3]

Product data:

Product Type: Primary Antibodies

Clone Name: 1A2E3
Applications: IHC, WB

Recommend Dilution: WB: 1 ug/mL, ICC: 2.5 ug/mL

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Mouse monoclonal MD-2 antibody was raised against a recombinant protein corresponding

to amino acids 2 to 160 of human MD-2

Formulation: MD-2 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.

Concentration: 1ug/ul

Purification: MD-2 Monoclonal Antibody is immunoaffinity chromotography purified IgG.

Gene Name: lymphocyte antigen 96

Database Link: NP 001182726 Entrez Gene 17087 MouseEntrez Gene 448830 RatEntrez Gene 23643 Human

Background: MD-2 Monoclonal Antibody: MD-2 is a member of the Toll/interleukin-1 receptor (TIR) family,

a group of proteins that include the Toll-like receptors (TLRs). TLRs are signaling molecules that recognize different pathogen-associated molecular patterns (PAMPs) and serve as an important link between the innate and adaptive immune responses. TLR4, the major signaling receptor for lipopolysaccharide (LPS), requires the binding of MD-2 to its extracellular region for maximal response to LPS. The specificity of this response is

determined by the species of MD-2; e.g., human MD-2 transfected into mouse cells can cause mouse TLR4 to react to LPS analogs that are normally antagonistic to human but not mouse

TLR4.

Synonyms: ESOP-1; ly-96; MD-2; MD2


Protein Families: Secreted Protein


Protein Pathways: Pathogenic Escherichia coli infection, Toll-like receptor signaling pathway





Product images:





Western blot analysis of (A) 25 and (B) 125 ng of MD-2 recombinant protein with MD-2 antibody at 1 ug/mL.

Immunohistochemistry of MD-2 in human spleen with MD-2 antibody at 2.5 ug/mL.