

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA324336**

MAP2 Mouse Monoclonal Antibody [Clone ID: 159CT34.12.3.4]

## **Product data:**

**Product Type:** Primary Antibodies

**Clone Name:** 159CT34.12.3.4

**Applications:** IF, IHC, WB

**Recommend Dilution:** WB: 1:1000, IF: 1:10~50, IHC: 1:10~50

Reactivity: Human Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: This MAP2 Monoclonal antibody was raised using purified His-tagged recombinant human

MAP2.

**Formulation:** Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.

**Concentration:** 0.35 mg/ml **Predicted Protein Size:** 199526 Da

Gene Name: microtubule associated protein 2

Database Link: NP 114034 Entrez Gene 4133 Human


Synonyms: DKFZp686I2148; MAP2A; MAP2B; MAP2C; microtubule-associated protein 2;

OTTHUMP00000163916


Protein Families: Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS



## **Product images:**




MAP2 Antibody (Cat. #TA324336) western blot analysis in MCF-7 cell line lysates (15ug/lane). This demonstrates the MAP2 antibody detected the MAP2 protein (arrow) (1:1000).



MAP2 Antibody (Ascites) (Cat. #TA324336)immunohistochemistry analysis in formalin fixed and paraffin embedded human brain tissue followed by peroxidase conjugation of the secondary antibody and DAB staining. This data demonstrates the use of MAP2 Antibody (Ascites) for immunohistochemistry. Clinical relevance has not been evaluated.





Confocal immunofluorescent analysis of MAP2 Antibody (Cat#TA324336) with brain tissue followed by Alexa Fluor 488-conjugated goat antimouse IgG (green). DAPI was used to stain the cell nuclear (blue).