

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA324352

Fatty Acid Synthase (FASN) Mouse Monoclonal Antibody [Clone ID: 497CT15.2.5]

Product data:

Product Type: Primary Antibodies

Clone Name: 497CT15.2.5

Applications: IF, WB

Recommend Dilution: IF: 1:200, WB: 1:100

Reactivity: Human, Mouse

Host: Mouse Isotype: IgM

Clonality: Monoclonal

Immunogen: This FASN antibody is generated from mice immunized with a KLH conjugated synthetic

peptide between 942-973 amino acids from the Central region of human FASN.

Formulation: Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.

Concentration: 2.0 mg/ml

Purification: This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.

Predicted Protein Size: 273427 Da

Gene Name: fatty acid synthase

Database Link: NP 004095 Entrez Gene 14104 MouseEntrez Gene 2194 Human


Synonyms: FAS; OA-519; SDR27X1 **Protein Families:** Druggable Genome

Protein Pathways: Fatty acid biosynthesis, Insulin signaling pathway, Metabolic pathways






Product images:



FASN Antibody (Center) (Cat. #TA324352) western blot analysis in mouse brain tissue lysates (35ug/lane). This demonstrates the FASN (Center) antibody detected the FASN (Center) protein (arrow).



IF image of HepG2 cells stained with FASN (Center) antibody. HepG2 cells were incubated with TA324352 FASN primary antibody (1:200, 2 h at RT). For secondary antibody, Alexa Fluor® 488 conjugated donkey anti-mouse antibody (green) was used (1:1000, 1h). Nuclei were counterstained with Hoechst 33342 (blue). Note the highly specific localization of the FASN immunosignal to the cytoplasm, supported by Human Protein Atlas Data (http://www.proteinatlas.org/ENSG00000169710).