

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA324501

TATA binding protein (TBP) Mouse Monoclonal Antibody [Clone ID: 830CT4.3.3]

Product data:

Product Type: Primary Antibodies

Clone Name: 830CT4.3.3

Applications: IF, WB

Recommend Dilution: WB 1:1000, IF 1:100

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Purified His-tagged TBP protein was used to produced this monoclonal antibody.

Formulation: Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.

Concentration: 0.37 mg/ml

Purification: This antibody is purified through a protein G column, followed by dialysis against PBS.

Predicted Protein Size: 37698 Da

Gene Name: TATA-box binding protein

Database Link: NP 003185 Entrez Gene 21374 MouseEntrez Gene 6908 Human

Synonyms: GTF2D; GTF2D1; HDL4; SCA17; TFIID


Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Basal transcription factors, Huntington's disease






Product images:




TBP Antibody (Cat. #TA324501) western blot analysis in Hela, HepG2, mouse NIH/3T3 cell line lysates (35ug/lane). This demonstrates the TBP antibody detected the TBP protein (arrow).



Western blot analysis of extracts (10ug) from 3 different cell lines by using anti-TBP monoclonal antibody (1:1000).





Immunofluorescent staining of HeLa cells using anti-TBP mouse monoclonal antibody (TA324501) (1:100).