

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA324585

# CDX1 Mouse Monoclonal Antibody [Clone ID: 937CT11.6.2]


## **Product data:**

Product Type:	Primary Antibodies
Clone Name:	937CT11.6.2
Applications:	WB
<b>Recommend Dilution:</b>	WB: 1:1000
Reactivity:	Human, Mouse
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Purified His-tagged CDX1 protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	0.3 mg/ml
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Predicted Protein Size:	28138 Da
Gene Name:	caudal type homeobox 1
Database Link:	<u>NP_001795 Entrez Gene 12590 MouseEntrez Gene 1044 Human</u>
Synonyms:	MGC116915
Protein Families:	Transcription Factors



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

### **Product images:**



CDX1 Antibody (C-term) (Cat. #TA324585) western blot analysis in WiDr, Hela cell line and mouse liver tissue lysates (35ug/lane).This demonstrates the CDX1 antibody detected the CDX1 protein (arrow).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US