

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA324831

USP22 Mouse Monoclonal Antibody [Clone ID: 1154CT13.1.3.1]

Product data:


Product Type:	Primary Antibodies
Clone Name:	1154CT13.1.3.1
Applications:	WB
Recommend Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Purified His-tagged USP22 protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	0.5 mg/ml
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Predicted Protein Size:	59961 Da
Gene Name:	ubiquitin specific peptidase 22
Database Link:	<u>NP_056091 Entrez Gene 23326 Human</u>
Synonyms:	USP3L
Protein Families:	Protease



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:



USP22 Antibody (C-term) (Cat. #TA324831) western blot analysis in Daudi cell line lysates (35ug/lane).This demonstrates the USP22 antibody detected the USP22 protein (arrow).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US