

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA324850

NLK Mouse Monoclonal Antibody [Clone ID: 1146CT24.2.1]

Product data:

Product Type: Primary Antibodies

Clone Name: 1146CT24.2.1

Applications: WB

Recommend Dilution: WB: 1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Purified His-tagged NLK protein was used to produced this monoclonal antibody.

Formulation: Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.

Concentration: 0.5 mg/ml

Purification: This antibody is purified through a protein G column, followed by dialysis against PBS.

Predicted Protein Size: 58283 Da

Gene Name: nemo like kinase

Database Link: NP 057315 Entrez Gene 51701 Human


Synonyms: DKFZp761G1211; FLJ21033

Protein Families: Druggable Genome, Protein Kinase, Transcription Factors

Protein Pathways: Adherens junction, MAPK signaling pathway, Wnt signaling pathway



Product images:



NLK Antibody (Cat. #TA324850) western blot analysis in SiHa, Hela cell line lysates (35ug/lane). This demonstrates the NLK antibody detected the NLK protein (arrow).