


Product datasheet for **TA324910**

GOLPH3 Mouse Monoclonal Antibody [Clone ID: 905CT9.1.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	905CT9.1.1
Applications:	WB
Recommend Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified His-tagged GOLPH3 protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	0.41 mg/ml
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Predicted Protein Size:	33811 Da
Gene Name:	golgi phosphoprotein 3
Database Link:	NP_071413 Entrez Gene 64083 Human
Synonyms:	GOPP1; GPP34; MIDAS; Vps74

[View online »](#)

Product images:


GOLPH3 Antibody (C-term) (Cat. #TA324910) western blot analysis in 293, MDA-MB-231 cell line lysates (35ug/lane). This demonstrates the GOLPH3 antibody detected the GOLPH3 protein (arrow).