

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA327732

Papillary renal cell carcinoma (PRCC) Mouse Monoclonal Antibody [Clone ID: PN-15]

Product data:

Product Type: Primary Antibodies

Clone Name: PN-15
Applications: IHC

Recommend Dilution: IHC: 1:50 - 1:200

Reactivity: Human Host: Mouse

Isotype: IgG1, kappa **Clonality:** Monoclonal

Formulation: The antibody is diluted in tris buffered saline, pH 7.3-7.7, with 1% BSA and <1% sodium azide.

Purification: Affinity purification

Gene Name: papillary renal cell carcinoma (translocation-associated)

Database Link: NP 005964 Entrez Gene 5546 Human

Synonyms: RCCP1; TPRC

Note: Anti-renal cell carcinoma (RCC) recognizes a 200 kD glycoprotein localized in the brush border


of the proximal renal tubule. This antibody immunoreacts with approximately 90% of primary renal cell carcinomas and approximately 85% of metastatic renal cell carcinomas. Therefore, anti-RCC is a reliable tool for differentiating primary or metastatic renal cell carcinoma from non-renal tumors. It may be utilized as a marker for the differential diagnosis of eosinophilic renal tumors-granular variant of renal cell carcinoma, chromophobe renal cell carcinoma, and oncocytoma. Other tumors that may react with this antibody are parathyroid adenoma and an occasional breast carcinoma. Nephroblastoma, oncocytoma, mesoblastic nephroma, transitional cell carcinoma, and angiomyolipoma are not labeled with this antibody.

Protein Families: Druggable Genome





Product images:



Immunohistochemistry staining of Paraffin Renal cell carcinoma tissue by Renal Cell Carcinoma antibody (dilution: 1:50 - 1:200; visualization of staining: Cytoplasmic, membranous)