

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA328045

H3FT (HIST3H3) Mouse Monoclonal Antibody [Clone ID: 11D8]

Product data:

Product Type: Primary Antibodies

Clone Name: 11D8
Applications: IF, WB
Recommend Dilution: WB, IF
Reactivity: Human
Host: Mouse

Isotype:IgG2b, kappaClonality:Monoclonal

Immunogen: Modified synthetic peptide conjugated to KLH

Formulation: This antibody is provided in phosphate-buffered solution, pH 7.2, containing 0.05% sodium

azide.

Concentration: 0.5 mg/ml

Purification: The antibody was purified by antigen-affinity chromatography.

Gene Name: histone cluster 3, H3

Database Link: NP 003484 Entrez Gene 8290 Human

Background: Histone H3 is phosphorylated at serine 10 during mitosis and found to be involved in

transcriptional activation, chromatin decondensation, and chromosome compaction during

cell division, by the action of Aurora kinase and NIMA kinases.


Synonyms: g; H3; H3.4; H3FT; H3t

Protein Pathways: Systemic lupus erythematosus






Product images:



293T cell extracts were resolved by electrophoresis, transferred to nitrocellulose, and probed with purified anti-H3-pS10 clone 11D8 antibody. Proteins were visualized using an antimouse IgG secondary conjugated to HRP and chemiluminescence detection. Lane 1, is inter phase 293T cell extract; lane 2, is M phase 293T cell extract (overnight nocodazole-treated cells).



Exponentially growing Hela cells were stained with purified phospor-Histone H3 (Ser10), clone 11D8 antibody, followed by Alexa Fluor® 488 conjugated to anti-mouse IgG and DAPI. The phosphor-Histone staining is specific on metaphase cells only.