

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA328117

Topoisomerase I (TOP1) Mouse Monoclonal Antibody [Clone ID: 1291CT875.142.166]

Product data:

Product Type: Primary Antibodies
Clone Name: 1291CT875.142.166

Applications: WB

Reactivity: WB: 1:1000

Human

Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: This TOP1 antibody is generated from a mouse immunized with a recombinant protein from

the N-terminal region of human TOP1.

Formulation: PBS with 0.09% (W/V) sodium azide

Concentration: 0.5mg/ml

Purification: This antibody is purified through a protein G column, followed by dialysis against PBS.

Gene Name: topoisomerase (DNA) I

Database Link: NP 003277 Entrez Gene 7150 Human

Background: Releases the supercoiling and torsional tension of DNA introduced during the DNA

replication and transcription by transiently cleaving and rejoining one strand of the DNA duplex. Introduces a single-strand break via transesterification at a target site in duplex DNA. The scissile phosphodiester is attacked by the catalytic tyrosine of the enzyme, resulting in the formation of a DNA-(3'-phosphotyrosyl)-enzyme intermediate and the expulsion of a 5'-OH DNA strand. The free DNA strand then undergoes passage around the unbroken strand thus removing DNA supercoils. Finally, in the religation step, the DNA 5'-OH attacks the covalent intermediate to expel the active-site tyrosine and restore the DNA phosphodiester backbone (By similarity). Regulates the alternative splicing of tissue factor (F3) pre-mRNA in endothelial

cells.


Synonyms: TOPI

Protein Families: Druggable Genome





Product images:



Western blot analysis of lysates from HUVEC, Jurkat, MCF-7, PC-12 cell line (from left to right), using TOP1 Antibody (N-term) (Cat. #TA328117). TA328117 was diluted at 1:1000 at each lane. A goat anti-mouse IgG H&L (HRP) at 1:5000 dilution was used as the secondary antibody. Lysates at 35ug per lane.