

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336420

NIFK Mouse Monoclonal Antibody [Clone ID: 18E148]

Product data:

Product Type: Primary Antibodies

Clone Name: 18E148

Applications: WB

Recommend Dilution: WB: 0.5-3 ug/ml

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: This antibody was developed against a synthetic peptide corresponding to a portion of

human NIFK between amino acids 100-200.

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Gene Name: nucleolar protein interacting with the FHA domain of MKI67

Database Link: NP 115766 Entrez Gene 84365 Human

Background: NIFK is a nucleolar protein that interacts with Ki67, a protein strongly associated with, and a

specific marker of, cell proliferation. The marsupial counterpart of human NIFK was recently

identified and suggests that NIFK may be involved in the organization of higher order

chromatin structure. NIFK was found to interact with the forkhead associated domain of Ki67

in mitotic cells and the interaction is thought to be phosphorylation-dependent.


Synonyms: MKI67IP; Nopp34

Protein Families: Transcription Factors





Product images:



Western Blot: NIFK Antibody (18E148) TA336420 - Analysis of NIFK in 30 ugs of Ramos cell lysate using NIFK antibody at 0.5 ug/ml.