

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336538

Apolipoprotein E (APOE) Mouse Monoclonal Antibody [Clone ID: WUE-4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	WUE-4
Applications:	FC, WB
Recommend Dilution:	WB: 2 ug/ml, ELISA: 1:100-1:2000, FC: 1 ug per million cells, IHC: 1:50-1:200, IP: 1:10-1:500
Reactivity:	Human, Mouse (Does not react with: Rat)
Host:	Mouse
lsotype:	lgG1, kappa
Clonality:	Monoclonal
Immunogen:	Purified human ApoE [UniProt# P02649]
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at - 20C long term. Avoid freeze-thaw cycles.
Concentration:	0.85 mg/ml
Purification:	Protein G purified
Predicted Protein Size:	36 kDa
Gene Name:	apolipoprotein E
Database Link:	<u>NP_000032 Entrez Gene 11816 MouseEntrez Gene 25728 RatEntrez Gene 348 Human</u>




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	polipoprotein E (APOE) Mouse Monoclonal Antibody [Clone ID: WUE-4] – TA336538
Background:	ApoE (apolipoprotein E) is the protein constituent of cholesterol/triglyceride-rich plasma lipoproteins, and is a multifunctional glycosylated secretary protein found almost in all organs with high activity in hepatic tissues. ApoE expression is induced by cholesterol-rich diets and is enhanced in lipoproteins in humans with genetic disorder type III hyperlipoproteinemia (HLP) characterized by remnant lipoproteins accumulation in plasma and premature atherosclerosis. ApoE circulates in blood as a protein component of VLDLs, chylomicron remnants, a subclass of HDL etc, and in cerebrospinal fluid as well as CNS interstitial fluid on small particles and disks resembling HDLs. ApoE facilitates transport of cholesterol and other lipids, as well as the clearance of plasma lipoproteins by serving as a critical ligand for lipoprotein uptake by LDL receptors and related proteins. ApoE participates in lipids redistribution to cells (e.g. CNS) that require cholesterol and phospholipids for reparative processes. ApoE also involves in proliferation inhibition of smooth muscle cells/lymphocytes, antigen presentation, and cholesterol efflux stimulation from foam cell macrophages. Defects in APOE have been linked to HLPP3 (hyperlipoproteinemia type 3), AD2 (Alzheimer disease type 2), SBHD (sea-blue histiocyte disease), LPG (lipoprotein glomerulopathy), and certain autoimmune disorders including multiple sclerosis and psoriasis.
Synonyms:	AD2; APO-E; LDLCQ5; LPG
Note:	This ApoE antibody is useful for Western blot, ELISA, Immunohistochemistry and Immunoprecipitation. In Western blot a band is observed at ~36 kDa, representing the ApoE protein.
Protein Families:	Adult stem cells, Druggable Genome, Secreted Protein, Stem cell - Pluripotency
Protein Pathways:	Alzheimer's disease
Product image	S:

ApoE antibody TEST Marker human mouse human mouse rat 150 kDa 100 kDa 75 kDa 50 kDa 37 kDa . 25 kDa

Western Blot: ApoE Antibody (WUE-4) [TA336538] - 50 ug protein per lane for liver lysate and 0.5 ul of plasma at 9% SDS. Samples are loaded onto the gel as the following order: human liver lysate, mouse liver lysate, rat liver lysate, human plasma, and mo

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US