

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA336584**

## c-Myc (MYC) Mouse Monoclonal Antibody [Clone ID: 9E11]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: 9E11
Applications: FC, WB

**Recommend Dilution:** WB: 1:500-1:1000, ChIP: 2 ug / 500 ug extract, ELISA: 1:100-1:2000, FC: 1:200-1:400, IHC: 1:100,

IHC-F: 1:100, IHC-P: 1:100, IP: 2ug/mg lysate

**Reactivity:** Chicken, Human, Mouse, Yeast

**Host:** Mouse

Isotype: IgG2a, kappa
Clonality: Monoclonal

Immunogen: A synthetic peptide corresponding to amino acids 408-420 (AEEQKLISEEDL) of human c-Myc,

conjugated to KLH. [UniProt# P01106]

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 1.1 mg/ml

**Purification:** Protein A purified

Gene Name: v-myc avian myelocytomatosis viral oncogene homolog

Database Link: NP 002458 Entrez Gene 17869 MouseEntrez Gene 4609 Human



Background:

Myc genes are a family of proto-oncogenes (L- Myc, N- Myc and C- Myc) that codes for Myc proteins which are transcriptor factors implicated in cellular proliferation, differentiation, apoptosis, metabolism, adhesion and self-renovation of tumor stem cells. Myc protein can act as transcriptional activator/repressor, and is activated via response to diverse mitogenic signals (including Wnt, Shh and EGF) and has been found to be up-regulated in several types of cancers. c-Myc participates gene transcription regulation and binds DNA in a non-specific manner, yet can specifically recognizes core sequence 5'-CAC[GA]TG-3' also. c-Myc heterodimerization with another bHLH protein namely Myc-associated factor X (MAX) is required for efficient c-Myc- DNA binding. c-Myc interacts with several proteins such as TAF1C, SPAG9, PARP10, KDM5A, KDM5B, NO66, PIM2 and with FBXW7 when phosphorylated at Thr-58/Ser-62. c-Myc activate the transcription of growth-related genes and c- Myc overexpression induce cell-cycle progression thereby implicating in a variety of cancers. Moreover, a chromosomal aberration involving c-Myc has been linked to a form of B-cell chronic lymphocytic leukemia and defective c-MYC is responsible for Burkitt lymphoma also.

**Synonyms:** bHLHe39; c-Myc; MRTL; MYCC

Note: This c-Myc antibody (clone 9E11) is useful for Flow Cytometry, ChIP, Immunoprecipitation,

ELISA, Immunohistochemistry- Frozen and Paraffin and Western blot.

**Protein Families:** Druggable Genome, Embryonic stem cells, Induced pluripotent stem cells, Stem cell -

Pluripotency, Stem cell relevant signaling - JAK/STAT signaling pathway, Stem cell relevant signaling - TGFb/BMP signaling pathway, Stem cell relevant signaling - Wnt Signaling pathway,


**Transcription Factors** 

Protein Pathways: Acute myeloid leukemia, Bladder cancer, Cell cycle, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Jak-STAT signaling pathway, MAPK

signaling pathway, Pathways in cancer, Small cell lung cancer, TGF-beta signaling pathway,


Thyroid cancer, Wnt signaling pathway

## **Product images:**



Western Blot: c-Myc Antibody (9E11) TA336584 - Analysis of c-Myc in Jurkat whole cell lysate.





Flow Cytometry: c-Myc Antibody (9E11) TA336584 - c-Myc antibody was tested at 1:400 in HL-60 cells using an Alexa Fluor 488 secondary (shown in purple). M1 is defined by unstained cells.