

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336637

IGF2R Mouse Monoclonal Antibody [Clone ID: 2G11]

Product data:

Product Type: Primary Antibodies

Clone Name: 2G11

Applications: FC, IF, IHC

Recommend Dilution: WB: 1 ug/ml, ELISA: , FC: 1 ug per million cells, IF: 1:10-1:250, IHC: 1:100, IHC-P: 1:100, IP: 1:10-

1:500

Reactivity: Bovine, Human, Primate, Rat

Host: Mouse

Isotype: IgG2a, kappa
Clonality: Monoclonal

Immunogen: Purified bovine Mannose 6 Phosphate Receptor (Cation independent) [UniProt# P08169]

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 1 mg/ml

Purification: Protein A purified

Gene Name: insulin like growth factor 2 receptor

Database Link: NP 000867 Entrez Gene 25151 RatEntrez Gene 3482 Human

Background: Mannose 6 phosphate is a sequence tag that plays a pivotal role in delivering proteins to the

lysosome. Specifically, the mannose 6 phosphate sequence is tagged to N-linked

oligosaccharides of lysosomal hydrolases as they pass through the cis-Golgi network. Once the lysosomal enzymes are tagged with the M6P sequence, they are destined to the late endosome by way of vesicular transport. After a low pH dissociates the M6P Receptor from its ligand, the M6P sequences are recycled. Additionally, M6P is bound by lectin in the immune system, and is converted to fructose 6-phosphate by mannose phosphate

isomerase.

Synonyms: CD222; CIMPR; M6P-R; MPR1; MPRI





Note: By Western Blot, a 300 kDa protein representing Mannose 6 Phosphate Receptor is seen in

HeLa cell extract under non-reducing conditions. Immunofluorescence staining of Mannose 6 Phosphate Receptor in HeLa cells with this antibody results in perinuclear staining. This


antibody is also useful in Flow Cytometry, Immunoprecipitation, ELISA and

Immunohistochemistry-on parraffin embedded sections.


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Lysosome

Product images:




Immunohistochemistry: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) TA336637 - Staining of human testis using TA336637 antibody at 1:100.




Immunocytochemistry/Immunofluorescence: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) TA336637 - HMVEC Cells






Immunocytochemistry/Immunofluorescence: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) TA336637 - Mannose 6 Phosphate Receptor (Cation independent) antibody (2G11) was tested at (1:250) in HeLa cells with Dylight 488 (green). Nuclei a




Immunocytochemistry/Immunofluorescence:
Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) TA336637 Mannose 6-Phosphate Receptor staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue) is shown. Cells were g



Immunocytochemistry/Immunofluorescence: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) TA336637 -Staining in Hela Cells. Mannose 6-Phosphate Receptor staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue





Flow Cytometry: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) TA336637 - Intracellular flow cytometric staining of 1 x 10^6 MCF-7 cells using Mannose 6 Phosphate Receptor antibody (dark blue). Isotype control shown in orange. An anti