

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA336671

### TRAF5 Mouse Monoclonal Antibody [Clone ID: 55A219]

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	55A219
Applications:	WB
<b>Recommend Dilution:</b>	WB: 2 ug/ml
Reactivity:	Human
Host:	Mouse
lsotype:	IgG
Clonality:	Monoclonal
Immunogen:	Anti-TRAF5 monoclonal antibody was raised against a fusion protein corresponding to amino acids 77 to 186 of human TRAF5.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at - 20C long term. Avoid freeze-thaw cycles.
Concentration:	0.5 mg/ml
Purification:	Protein G purified
Gene Name:	TNF receptor associated factor 5
Database Link:	<u>NP_665702 Entrez Gene 7188 Human</u>
Background:	Tumor necrosis factor (TNF) induced signaling is mediated through association of TNF receptor (TNFR) with adaptor proteins, such as TNF receptor associated proteins (TRAFs). TRAFs form a family of cytoplasmic adapter proteins that mediate signal transduction from many members of the TNF-receptor superfamily (e.g. RANK, CD30, CD40, etc.) and the interleukin-1 receptor. The carboxy-terminal region of TRAFs is required for self-association and interaction with receptor cytoplasmic domains following ligand-induced oligomerization. Recent molecular cloning studies have lead to identification of six TRAFs (TRAF1-TRAF6). TRAF5 is a 558-amino acid protein. TRAF5 is implicated in NF-kB and c-Jun NH(2)-terminal kinase/stress-activated protein kinase activation by members of the TNF receptor superfamily, including CD27, CD30, CD40, and lymphotoxin-beta receptor . Targeted disruption of TRAF5 gene causes defects in CD40-CD27 mediated lymhocyte activation.
Synonyms:	MGC:39780; RNF84




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	TRAF5 Mouse Monoclonal Antibody [Clone ID: 55A219] – TA336671
Note:	This antibody can be used in Western Blot where a 66 kDa band is seen (predicted molecular weight 66.8 kDa).
Protein Families	: Druggable Genome

Protein Pathways: Pathways in cancer, Small cell lung cancer

## **Product images:**



Western Blot: TRAF5 Antibody (55A219) TA336671 - Analysis of TRAF5 using this antibody at 2 ug/ml dilution against 10 ug of HeLa cell lysate.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US