

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336673

Noxa (PMAIP1) Mouse Monoclonal Antibody [Clone ID: 114C307.1]

Product data:


Product Type:	Primary Antibodies
Clone Name:	114C307.1
Applications:	IF, IHC, WB
Recommend Dilution:	WB: 1-2 ug/ml, IF: 1:10-1:500, IHC: 1:10-1:500, IHC-P: 5 ug/ml
Reactivity:	Human, Mouse
Host:	Mouse
lsotype:	lgG1, kappa
Clonality:	Monoclonal
Immunogen:	Noxa Antibody was developed by immunizing mice with a fusion protein containing human Noxa.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at - 20C long term. Avoid freeze-thaw cycles.
Concentration:	0.5 mg/ml
Purification:	Protein G purified
Gene Name:	phorbol-12-myristate-13-acetate-induced protein 1
Database Link:	<u>NP_066950 Entrez Gene 58801 MouseEntrez Gene 5366 Human</u>
Background:	Recently, a BH-3 only member of the Bcl-2 family have been identified in human and mouse and designated as Noxa (for damage) (1). The expression of the Noxa gene involves direct activation of its promoter by p53. Increased expression of Noxa protein occurs in normal thymocytes but not in p53-deficient thymocytes. Noxa cDNA codes for a 103-amino acid protein. The coimmunoprecipitation data suggest that Noxa protein may interact with proteins belonging to the Bcl-2 family, such as, Bcl-XL and Mcl-1. Oda et al (1) have also shown that blocking the endogenous Noxa induction results in the suppression of apoptosis. Treatment of cells with Noxa antisense oligonucleotides blocks radiation-induced apoptosis suggesting that Noxa may represent a mediator of p53-dependent apoptosis.
Synonyms:	APR; NOXA
Note:	Immunocytochemistry and Immunofluorescence; See Yakovlev et all (2003) for details.
Protein Families:	Druggable Genome, Stem cell - Pluripotency




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US Protein Pathways:

p53 signaling pathway


Product images:



Western Blot: Noxa Antibody (114C307.1) [TA336673] - detection of Noxa antibody in A) 293 cells transfected with human Noxa, B) RL-7 cell (a follicular lymphoma) and C) T98G lysate using this antibody.



Immunohistochemistry-Paraffin: Noxa Antibody (114C307.1) [TA336673] - Bouin-fixed, paraffinembedded Chronic Lymphocytic Leukemia (CLL) xenograft stained with Noxa antibody (1:2000), peroxidase-conjugate and DAB chromogen. Staining for 2 hr at RT was used

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

Immunocytochemistry/Immunofluorescence: Noxa Antibody (114C307.1) [TA336673] - Noxa antibody was tested in HeLa cells with Dylight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and Dylight 550 (red). An antibody dilution of 1:

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US