

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336873

TLR5 Mouse Monoclonal Antibody [Clone ID: 85B152.5]

Product data:

Product Type: Primary Antibodies

Clone Name: 85B152.5

Applications: FC, WB

Recommend Dilution: WB: 1-3 ug/ml, FC: 1-3 ug/10^6 cells in 100 ul

Reactivity: Canine, Human, Mouse

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: This antibody was developed against KLH-conjugated synthetic peptide corresponding to a

portion of human TLR5. It will cross-react with mouse TLR5 (NP 003259).

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Predicted Protein Size: 100 kDa

Gene Name: toll like receptor 5

Database Link: NP 003259 Entrez Gene 53791 MouseEntrez Gene 7100 Human

Background: The Toll-like receptor (TLR) family in mammal comprises a family of transmembrane proteins

characterized by multiple copies of leucine rich repeats in the extracellular domain and IL-1 receptor motif in the cytoplasmic domain. Like its counterparts in Drosophila, TLRs signal through adaptor molecules and could constitute an important and unrecognized component of innate immunity in humans. The TLR family is a phylogenetically conserved mediator of innate immunity that is essential for microbial recognition. TLRs characterized so far activate the MyD88/interleukin-1 receptor-associated kinase (IRAK) signaling pathway. Thirteen homologs of TLRs (TLR1-13) have been described. Toll-like receptor 5 (TLR5) expression is upregulated following exposure to bacteria or to the TLR5 agonist, flagellin. Gram-negative bacteria, stimulate monocyte/macrophage cells in a TLR5-specific, CD14-independent manner. The TLR5 receptor thus appears to be the principal means by which the innate

immune system recognizes flagellated bacterial pathogens.



TLR5 Mouse Monoclonal Antibody [Clone ID: 85B152.5] - TA336873


Synonyms: MELIOS; SLE1; SLEB1; TIL3

Note: Reported use in Flow Cytometry (cell surface): see Wong et al, 2007 for details.


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Pathogenic Escherichia coli infection, Toll-like receptor signaling pathway

Product images:



Western Blot: TLR5 Antibody (85B152.5) TA336873 - Normal human bronchial epithelial total cell lysate. Image from verified customer review.



Flow Cytometry: TLR5 Antibody (85B152.5) TA336873 - Intracellular flow analysis of TLR5 in Ramos cells using 1 ug of TLR5 antibody. Shaded histogram represents Ramos cells without antibody; green represents isotype control antibody; red represents TLR5