

# Product datasheet for TA337048

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Bif (SH3GLB1) Mouse Monoclonal Antibody [Clone ID: 30A882.1.1]

### **Product data:**

**Product Type:** Primary Antibodies

**Clone Name:** 30A882.1.1

Applications: WB

Recommend Dilution: WB: 1-2 ug/ml, IHC: 5ug/ml, IP: 1:10-1:500

Reactivity: Human Host: Mouse

Isotype:IgG1, kappaClonality:Monoclonal

**Immunogen:** A full-length human Bif-1 recombinant protein was used as immunogen.

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

**Concentration:** 0.5 mg/ml

**Purification:** Protein G purified

Gene Name: SH3 domain containing GRB2 like endophilin B1


Database Link: NP 057093 Entrez Gene 51100 Human
Synonyms: Bif-1; CGI-61; dJ612B15.2; PPP1R70

**Protein Pathways:** Endocytosis





## **Product images:**



Western Blot: Bif-1/Endophilin-B1 Antibody (30A882.1.1) TA337048 - Analysis for Bi-f-1/Endophilin-B1 using this antibody at 2 ug/ml dilution against 10 ug of Jurkat cell lysate.