

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA337074

p73 (TP73) Mouse Monoclonal Antibody [Clone ID: 38C674.2]

Product data:

Product Type: Primary Antibodies

Clone Name: 38C674.2 Applications: IHC, WB

Recommend Dilution: WB: 1-3 ug/ml, ChIP: 5ug, IF: 1:10-1:500, IHC-F: 1:10-1:500, IHC-P: 5 ug/ml

Reactivity: Human, Mouse, Rat

Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: This antibody was developed against a peptide corresponding to a portion of amino acid

residues 1-50 of human dNp73. Mouse and human sequences are more than 95 % identical

at these amino acid residues. In mouse, alanine at position 8 is replaced by methio

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Gene Name: tumor protein p73

Database Link: NP 001119712 Entrez Gene 362675 RatEntrez Gene 7161 Human

Synonyms: P73

Note: Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM sodium

citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 min.


Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Neurotrophin signaling pathway, p53 signaling pathway







Product images:




Western Blot: Delta Np73 Antibody (38C674.2) TA337074 - Delta Np73 Antibody (38C674.2) analysis for Np73 using TA337074 at 1 ug/ml in A) a cell line transfected with Np73 cDNA, B) HeLa, and C) NIH 3T3 cell lysate.






Immunohistochemistry-Frozen: Delta Np73 Antibody (38C674.2) TA337074 - Delta Np73 Antibody (38C674.2) analysis of Delta Np73 in embryonic day 15 (A) and day 17 (B) E15 rat embryonic kidney at 1:2000. Delta Np73 is expressed in ureteric bud (UB) branches





Immunohistochemistry-Paraffin: Delta Np73 Antibody (38C674.2) TA337074 - Delta Np73 Antibody (38C674.2) Formalin-fixed, paraffinembedded lung squamous cell carcinoma probed with p73 antibody at 5 ug/ml.



Immunohistochemistry-Paraffin: Delta Np73 Antibody (38C674.2) TA337074 - Delta Np73 Antibody (38C674.2) Staining of human placenta. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentra