

Product datasheet for **TA337076**

KCNIP1 Mouse Monoclonal Antibody [Clone ID: 6C239.4]


Product data:

Product Type:	Primary Antibodies
Clone Name:	6C239.4
Applications:	WB
Recommend Dilution:	WB: 1-3 ug/ml
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Full length recombinant protein was used as the immunogen.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	0.5 mg/ml
Purification:	Protein G purified
Gene Name:	potassium voltage-gated channel interacting protein 1
Database Link:	NP_055407 Entrez Gene 30820 Human
Synonyms:	KCHIP1; VABP
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)

Product images:



Western Blot: VABP Antibody (6C239.4) TA337076
- Analysis of VABP in 1) transfected 293T cells from Abnova and 2) untransfected 293T cells using TA337076 at 2 ug/ml.