

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA337154

IL17E (IL25) Mouse Monoclonal Antibody [Clone ID: 8B1]

Product data:

Product Type: Primary Antibodies

Clone Name: 8B1

Applications: FC, WB

Recommend Dilution: WB: 0.5-2ug/ml~, FC: 1-2ug/ml

Reactivity: Human Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: Full-length human IL17E recombinant protein was used as the immunogen for this antibody.

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Gene Name: interleukin 25

Database Link: NP 758525 Entrez Gene 64806 Human

Background: As a member of the IL-17 family, IL-25 has been identified as a pro-inflammatory cytokine

which is produced by Th17 effector cells thus enabling the TH17 subset to initiate or maintain the autoimmune process in diseases such as Rheumatoid Arthritis (1-3). IL-25 appears to stimulate IL-8 production and itself appears to be activated through IL-23, IL-6, and TGF-Beta (1). Studies of IL-25 expression and activation through NF-kB taken together with IL-17 production and Th17 polarization through TLR ligand activation further closely tie the IL-17 and Th17 family together with the TLR signaling network initiating and controlling immune

responses (4-7).

Synonyms: IL17E


Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction






Product images:



Western Blot: IL17E Antibody (8B1) TA337154 - Analysis of IL17E in full-length recombinant protein using IL17E antibody at 0.5 ug/ml.



Flow Cytometry: IL17E Antibody (8B1) TA337154 - Intracellular analysis of IL-17E in PC3 cells using 1 ug of IL-17E antibody (green histogram represents isotype control antibody this antibody, red represents IL-17E antibody. this antibody was used for t