

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA346932

Pyruvate Dehydrogenase E2 (DLAT) Mouse Monoclonal Antibody [Clone ID: 4A4-B6-C10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	4A4-B6-C10
Applications:	IF, IP, WB
Recommend Dilution:	WB: 1:1000, IF: 1:300
Reactivity:	Human, Mouse
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for DLAT antibody: purified recombinant human Pyruvate Dehydrogenase E2 protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.02% sodium azide, 50%,glycerol
Purification:	Affinity purified
Predicted Protein Size:	69 kDa
Predicted Protein Size: Gene Name:	69 kDa dihydrolipoamide S-acetyltransferase
Gene Name:	dihydrolipoamide S-acetyltransferase




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Pyruvate Dehydrogenase E2 (DLAT) Mouse Monoclonal Antibody [Clone ID: 4A4-B6-C10] –
	TA346932

Protein Families: Druggable Genome

Protein Pathways: Citrate cycle (TCA cycle), Glycolysis / Gluconeogenesis, Metabolic pathways, Pyruvate metabolism

Product images:




Western blot detection of Pyruvate Dehydrogenase E2 in Jurkat, A549, U251, F9, Lncap and Hela cell lysates using Pyruvate Dehydrogenase E2 mouse mAb (1:1000 diluted).Predicted band size:69KDa.Observed band size:69KDa.

Immunocytochemistry stain of Hela using Pyruvate Dehydrogenase E2 mouse mAb (1:300).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Pyruvate Dehydrogenase E2 (DLAT) Mouse Monoclonal Antibody [Clone ID: 4A4-B6-C10] – TA346932



Immunoprecipitation analysis of Hela cell lysates using Pyruvate Dehydrogenase E2 mouse mAb.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US