

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA347009

LCN2 Mouse Monoclonal Antibody [Clone ID: 4C12-C9-A10]

Product data:

Product Type: Primary Antibodies

Clone Name: 4C12-C9-A10

Applications:ELISAReactivity:HumanHost:MouseIsotype:IgG2b

Clonality: Monoclonal

Immunogen: The immunogen for NGAL antibody: Recombinant protein of human NGAL.

Formulation: Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and

50% glycerol.

Purification: Affinity purified

Gene Name: lipocalin 2

Database Link: NP 005555 Entrez Gene 3934 Human

Background: Lipocalin-2 (LCN2), also known as neutrophil gelatinase-associated lipocalin (NGAL), is a

protein that in humans is encoded by the LCN2 gene. NGAL is involved in innate immunity by sequestrating iron that in turn limits bacterial growth. It is expressed in neutrophils and in low levels in the kidney, prostate, and epithelia of the respiratory and alimentary tracts. In the case of acute kidney injury (AKI), NGAL is secreted in high levels into the blood and urine


within 2 hours of injury. So NGAL is used as a biomarker of kidney injury.

Synonyms: 24p3; MSFI; NGAL **Protein Families:** Secreted Protein





Product images:



Indirect ELISA assay for anti-NGAL mouse mAb.Antigen coating concentration: 4ug/ml.