

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA347030

PRMT1 Mouse Monoclonal Antibody [Clone ID: 4E9-H9-H12]

Product data:

Product Type: Primary Antibodies

Clone Name: 4E9-H9-H12

Applications: WB

Recommend Dilution: WB: 1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: The immunogen for PRMT1 antibody: purified recombinant human PRMT1 protein fragments

expressed in E.coli.

Formulation: Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and

50% glycerol.

Purification: Affinity purified

Predicted Protein Size: 42 kDa

Gene Name: protein arginine methyltransferase 1

Database Link: NP 001527 Entrez Gene 3276 Human

Background: This gene encodes a member of the protein arginine N-methyltransferase (PRMT) family.Post-

translational modification of target proteins by PRMTs plays an important regulatory role in many biological processes, whereby PRMTs methylate arginine residues by transferring methyl groups from S-adenosyl-L-methionine to terminal guanidino nitrogen atoms. The encoded protein is a type I PRMT and is responsible for the majority of cellular arginine methylation activity. Increased expression of this gene may play a role in many types of cancer. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of


chromosome 5. [provided by RefSeq,Dec 2011]

Synonyms: ANM1; HCP1; HRMT1L2; IR1B4





Product images:



Western blot detection of PRMT1 in Hela, A549 and SW480 cell lysates using PRMT1 mouse mAb (1:1000 diluted).Predicted band size:42KDa.Observed band size:42KDa.