

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA347053

SMYD1 Mouse Monoclonal Antibody [Clone ID: 2E11-B6-C7]

Product data:

Product Type: Primary Antibodies

Clone Name: 2E11-B6-C7

Applications: WB

Recommend Dilution: WB: 1:100

Reactivity: Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: The immunogen for SMYD1 antibody: purified recombinant human SMYD1 protein fragments

expressed in E.coli.

Formulation: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with

0.02% sodium azide, 50% glycerol

Purification: Affinity purified

Predicted Protein Size: 56 kDa

Gene Name: SET and MYND domain containing 1

Database Link: NP 938015 Entrez Gene 297333 Rat

Background: Acts as a transcriptional repressor. Essential for cardiomyocyte differentiation and cardiac


morphogenesis By similarity.

Synonyms: BOP; KMT3D; ZMYND18; ZMYND22

Protein Families: Druggable Genome



Product images:



Western blot detection of SMYD1 in rat skeletal muscle lysates and using SMYD1 mouse mAb (1:100 diluted). Predicted band size: 56KDa.Observed band size: 56KDa.