

## Product datasheet for **TA347055**

### **SMYD4 Mouse Monoclonal Antibody [Clone ID: 7C8-B11-F10]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	7C8-B11-F10
Applications:	WB
Recommend Dilution:	WB: 1:200
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for SMYD4 antibody: purified recombinant human SMYD4 protein fragments expressed in E.coli.
Formulation:	ascites
Predicted Protein Size:	89 kDa
Gene Name:	SET and MYND domain containing 4
Database Link:	<a href="#">NP_443160</a> <a href="#">Entrez Gene 114826</a> <a href="#">Human</a>
Synonyms:	ZMYND21




[View online »](#)

## Product images:



Western blot detection of SMYD4 in HeLa, MCF7 and A431 cell lysates using SMYD4 mouse mAb (1:200 diluted). Predicted band size: 89KDa. Observed band size: 89KDa.



Western blot detection of SMYD4 in 40ug HeLa, 40ug MCF7 and 40ug MDA-MB-468 cell lysates using SMYD4 mouse mAb (1:100 diluted). Predicted band size: 89KDa. Observed band size: 89KDa.