

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500001

AKT3 Mouse Monoclonal Antibody [Clone ID: OTI9B2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI9B2

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:500, IHC 1:50, IF 1:100, FLOW 1:100

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG3

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human AKT3 (NP_005456) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.25 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 55.6 kDa

Gene Name: AKT serine/threonine kinase 3

Database Link: NP 005456 Entrez Gene 23797 MouseEntrez Gene 10000 Human

Background: The protein encoded by this gene is a member of the AKT, also called PKB, serine/threonine

protein kinase family. AKT kinases are known to be regulators of cell signaling in response to

insulin and growth factors. They are involved in a wide variety of biological processes including cell proliferation, differentiation, apoptosis, tumorigenesis, as well as glycogen synthesis and glucose uptake. This kinase has been shown to be stimulated by platelet-derived growth factor (PDGF), insulin, and insulin-like growth factor 1 (IGF1). Alternatively

splice transcript variants encoding distinct isoforms have been described.

Synonyms: MPPH; MPPH2; PKB-GAMMA; PKBG; PRKBG; RAC-gamma; RAC-PK-gamma; STK-2

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase







Protein Pathways:

Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Tight junction, Toll-like receptor signaling pathway, VEGF signaling pathway

Product images:







HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AKT3 ([RC221051], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AKT3. Positive lysates [LY401675] (100ug) and [LC401675] (20ug) can be purchased separately from OriGene.


AKT3 antibody (9B2) at 1:500 dilution + BALB/C 3T3 cell lysate





Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-AKT3 monoclonal antibody at 1:200 dilution.









Figure from citation: Western Blot of Akt3 protein level by using anti-Akt3 antibody in rat endometrium during pregnancy and pseudopregnancy. View Citation




Immunohistochemical staining of paraffinembedded Human lung tissue within the normal limits using anti-AKT3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500001)





Anti-AKT3 mouse monoclonal antibody (TA500001) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY AKT3 ([RC221051]).



HEK293T cells transfected with either [RC221051] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-AKT3 antibody (TA500001), and then analyzed by flow cytometry.