

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500074

Apolipoprotein M (APOM) Mouse Monoclonal Antibody [Clone ID: OTI3H3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3H3
Applications: IHC, WB

Recommend Dilution: WB: 1:1000, IHC 1;50

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Recombinant protein expressed in E.coli corresponding to amino acids 23-188 of human

APOM

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.85 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 21.1 kDa

Gene Name: apolipoprotein M

Database Link: NP 061974 Entrez Gene 55937 Human

Background: The APOM is an apolipoprotein and member of the lipocalin protein family. It is found

associated with high density lipoproteins and to a lesser extent with low density lipoproteins and triglyceride-rich lipoproteins. The encoded protein is secreted through the plasma membrane but remains membrane-bound, where it is involved in lipid transport. Two transcript variants encoding two different isoforms have been found for this gene, but only

one of them has been fully characterized


Synonyms: apo-M; G3a; HSPC336; NG20

Protein Families: Secreted Protein






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY APOM ([RC204679], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-APOM. Positive lysates [LY402735] (100ug) and [LC402735] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 1 cell line by using anti-APOM monoclonal antibody at 1:200 dilution.





Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-APOM monoclonal antibody (1:200).



Immunohistochemical staining of paraffinembedded liver within the normal limits using anti-APOM (TA500074) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).