

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA500400**

## Melanoma gp100 (PMEL) Mouse Monoclonal Antibody [Clone ID: OTI9G3]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI9G3

**Applications:** IF, IHC, WB

**Recommend Dilution:** WB 1:2000, IHC: 1:150

**Reactivity:** Human, Dog

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human SILV expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 1.9 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 70.3 kDa

**Gene Name:** premelanosome protein

Database Link: NP 008859 Entrez Gene 481102 DogEntrez Gene 6490 Human

**Background:** Could be a melanogenic enzyme. Could represent an oncofetal self-antigen that is normally

expressed at low levels in quiescent adult melanocytes but overexpressed by proliferating neonatal melanocytes and during tumor growth. Release of the soluble form, ME20-S, could

protect tumor cells from antibody mediated immunity


**Synonyms:** D12S53E; gp100; ME20; ME20-M; ME20M; P1; P100; PMEL17; SI; SIL; SILV

**Protein Families:** Secreted Protein, Transmembrane







### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SILV ([RC200663], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SILV. Positive lysates [LY402062] (100ug) and [LC402062] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-SILV monoclonal antibody.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-SILV mouse monoclonal antibody. (TA500400)






Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-SILV mouse monoclonal antibody. (TA500400)




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-SILV mouse monoclonal antibody. (TA500400)




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-SILV mouse monoclonal antibody. (TA500400)






Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-SILV mouse monoclonal antibody. (TA500400)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-SILV mouse monoclonal antibody. (TA500400)




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-SILV mouse monoclonal antibody. (TA500400)





Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-SILV mouse monoclonal antibody. (TA500400)



Anti-SILV mouse monoclonal antibody (TA500400) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY SILV ([RC200663]).