

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA500435**

## Rad9 (RAD9A) Mouse Monoclonal Antibody [Clone ID: OTI5D9]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5D9

Applications: IHC, WB

Recommend Dilution: WB 1:2000 IHC 1:50

Reactivity: Human, Monkey, Dog

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Full-length protein expressed in 293T cell transfected with human RAD9A expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1.2 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 42.5 kDa

Gene Name: RAD9 checkpoint clamp component A

Database Link: NP 004575 Entrez Gene 483696 DogEntrez Gene 712345 MonkeyEntrez Gene 5883 Human

**Background:** This gene product is highly similar to Schizosaccharomyces pombe rad9, a cell cycle

checkpoint protein required for cell cycle arrest and DNA damage repair in response to DNA damage. This protein is found to possess 3' to 5' exonuclease activity, which may contribute to its role in sensing and repairing DNA damage. It forms a checkpoint protein complex with RAD1 and HUS1. This complex is recruited by checkpoint protein RAD17 to the sites of DNA damage, which is thought to be important for triggering the checkpoint-signaling cascade.

Use of alternative polyA sites has been noted for this gene.


Synonyms: RAD9

**Protein Families:** Druggable Genome, Stem cell - Pluripotency







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RAD9A ([RC204439], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RAD9A. Positive lysates [LY401451] (100ug) and [LC401451] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-anti-RAD9Amonoclonal antibody.



Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-RAD9Amouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500435, Dilution 1:50)





Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-RAD9Amouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500435, Dilution 1:50)