

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA500465**

## IRF3 Mouse Monoclonal Antibody [Clone ID: OTI5D2]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5D2
Applications: IHC, WB

Recommend Dilution: WB 1:2000 IHC 1:50

Reactivity: Human, Monkey

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human IRF3 expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 47.2 kDa

**Gene Name:** interferon regulatory factor 3

Database Link: NP 001562 Entrez Gene 718970 MonkeyEntrez Gene 3661 Human

**Background:** IRF3 encodes interferon regulatory factor 3, a member of the interferon regulatory

transcription factor (IRF) family. IRF3 is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a complex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons alpha and beta, as well as other

interferon-induced genes.

Synonyms: interferon regulatory factor 3; MGC94729 **Protein Families:** Druggable Genome, Transcription Factors


Protein Pathways: Cytosolic DNA-sensing pathway, RIG-I-like receptor signaling pathway, Toll-like receptor

signaling pathway







### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IRF3 ([RC209951], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IRF3 (1:500).




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-anti-IRF3monoclonal antibody.




Western blot analysis of extracts (35ug) from 4 different cell lines and human plasma by using anti-IRF3 monoclonal antibody (1:500).






Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and IRF3-Knockout HeLa cells (KO, Cat# [LC833375]) were separated by SDS-PAGE and immunoblotted with anti-IRF3 monoclonal antibody TA500465 (1:500`). Then the blotted membrane was stripped and reprobed with anti-PCNA antibody as a loading control.




Immunohistochemical staining of paraffinembedded breast tissue within the normal limits using anti-IRF3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-IRF3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-IRF3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of kidney tissue using anti-IRF3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)




Immunohistochemical staining of paraffinembedded lung tissue within the normal limits using anti-IRF3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-IRF3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-IRF3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)




Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-IRF3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Carcinoma of pancreas tissue using anti-IRF3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)





Immunohistochemical staining of paraffinembedded thyroid tissue within the normal limits using anti-IRF3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-IRF3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)










Immunohistochemical staining of paraffinembedded Carcinoma of prostate tissue using anti-IRF3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-IRF3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500465, Dilution 1:50)