

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500505

Kindlin 2 (FERMT2) Mouse Monoclonal Antibody [Clone ID: OTI14A11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI14A11

Applications: IF, IHC, IP, WB

Recommend Dilution: WB 1:1000~2000, IHC 1:50, IF 1:50~100, IP 2ug/500ul

Reactivity: Human, Monkey, Rat, Dog

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human FERMT2 expression vector

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 77.9 kDa

Gene Name: fermitin family member 2

Database Link: NP 006823 Entrez Gene 289992 RatEntrez Gene 490693 DogEntrez Gene 693958

MonkeyEntrez Gene 10979 Human

Background: Participates in the connection between ECM adhesion sites and the actin cytoskeleton and

also in the orchestration of actin assembly and cell shape modulation. Recruits migfilin


(FBLP1) protein to cell-ECM focal adhesion sites

Synonyms: KIND2; mig-2; MIG2; PLEKHC1; UNC112; UNC112B






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FERMT2 ([RC205727], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FERMT2. Positive lysates [LY416377] (100ug) and [LC416377] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-FERMT2 monoclonal antibody.






Western blot analysis of extracts (10ug) from a mouse cell line and a mouse tissue by using anti-FERMT2 monoclonal antibody (1:200).




Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and FERMT2-Knockout HeLa cells (KO, Cat# [LC812552]) were separated by SDS-PAGE and immunoblotted with anti-FERMT2 monoclonal antibody TA500505 (1:500`). Then the blotted membrane was stripped and reprobed with anti-PCNA antibody as a loading control.




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-FERMT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500505)






Anti-FERMT2 mouse monoclonal antibody (TA500505) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FERMT2 ([RC205727]).



Immunofluorescent staining of HepG2 cells using anti-FERMT2 mouse monoclonal antibody (TA500505).



Immunoprecipitation (IP) of FERMT2 by using TrueMab monoclonal anti-FERMT2 antibodies (Negative control: IP without adding anti-FERMT2 antibody.). For each experiment, 500ul of DDK tagged FERMT2 overexpression lysates (at 1:5 dilution with HEK293T lysate), 2ug of anti-FERMT2 antibody and 20ul (0.1mg) of goat anti-mouse conjugated magnetic beads were mixed and incubated overnight. After extensive wash to remove any non-specific binding, the immunoprecipitated products were analyzed with rabbit anti-DDK polyclonal antibody.