

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA500533**

## **BubR1 (BUB1B) Mouse Monoclonal Antibody [Clone ID: OTI5D9]**

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI5D9

**Applications:** IF, IHC, WB

Recommend Dilution: WB 1:2000 IHC 1:50 IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human BUB1B expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1.2 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 119.5 kDa

**Gene Name:** BUB1 mitotic checkpoint serine/threonine kinase B

Database Link: NP 001202 Entrez Gene 701 Human

**Background:** This gene encodes a kinase involved in spindle checkpoint function. The protein has been

localized to the kinetochore and plays a role in the inhibition of the anaphase-promoting

complex/cyclosome(APC/C), delaying the onset of anaphase and ensuring proper

chromosome segregation. Impaired spindle checkpoint function has been found in many

forms of cancer.

Synonyms: Bub1A; BUB1beta; BUBR1; hBUBR1; MAD3L; MVA1; SSK1


**Protein Families:** Druggable Genome, Protein Kinase, Stem cell - Pluripotency

**Protein Pathways:** Cell cycle







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BUB1B ([RC205844], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BUB1B. Positive lysates [LY420068] (100ug) and [LC420068] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-BUB1B monoclonal antibody.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-BUB1B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-BUB1Bmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of kidney tissue using anti-BUB1Bmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of liver tissue using anti-BUB1Bmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Carcinoma of lung tissue using anti-BUB1Bmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-BUB1Bmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of pancreas tissue using anti-BUB1Bmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-BUB1Bmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of prostate tissue using anti-BUB1Bmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-BUB1Bmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500533, Dilution 1:50)





Anti-BUB1B mouse monoclonal antibody (TA500533) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY BUB1B ([RC205844]).