

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500554

Aurora C (AURKC) Mouse Monoclonal Antibody [Clone ID: OTI10A7]

Product data:


Primary Antibodies
OTI10A7
IF, WB
WB 1:500, IF 1:100
Human
Mouse
lgG1
Monoclonal
Full-length protein expressed in 293T cell transfected with human AURKC expression vector
PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
0.85 mg/ml
Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
35.6 kDa
aurora kinase C
<u>NP_001015878 Entrez Gene 6795 Human</u>
This gene encodes a member of the Aurora subfamily of serine/threonine protein kinases. The encoded protein is a chromosomal passenger protein that forms complexes with Aurora- B and inner centromere proteins and may play a role in organizing microtubules in relation to centrosome/spindle function during mitosis. This gene is overexpressed in several cancer cell lines, suggesting an involvement in oncogenic signal transduction. Alternative splicing results in multiple transcript variants.
AIE2; AIK3; ARK3; AurC; aurora-C; HEL-S-90; SPGF5; STK13
Druggable Genome, Protein Kinase




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:



Anti-AURKC mouse monoclonal antibody (TA500554) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY AURKC ([RC213890]).



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AURKC ([RC213890], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AURKC. Positive lysates [LY423127] (100ug) and [LC423127] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US