

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500650

XRCC4 Mouse Monoclonal Antibody [Clone ID: OTI4H9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4H9
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human XRCC4 expression vector

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 38.3 kDa

Gene Name: X-ray repair cross complementing 4

Database Link: NP 072044 Entrez Gene 7518 Human

Background: The protein encoded by this gene functions together with DNA ligase IV and the DNA-

dependent protein kinase in the repair of DNA double-strand break by non-homologous end joining and the completion of V(D)J recombination events. The non-homologous end-joining pathway is required both for normal development and for suppression of tumors. This gene functionally complements XR-1 Chinese hamster ovary cell mutant, which is impaired in DNA double-strand breaks produced by ionizing radiation and restriction enzymes. Alternative transcription initiation and alternative splicing generates several transcript variants.

Synonyms: SSMED


Protein Families: Druggable Genome

Protein Pathways: Non-homologous end-joining






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY XRCC4 ([RC218029], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-XRCC4. Positive lysates [LY402930] (100ug) and [LC402930] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-XRCC4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)




Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-XRCC4mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)




Immunohistochemical staining of paraffinembedded lung tissue within the normal limits using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-XRCC4mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)




Immunohistochemical staining of paraffinembedded thyroid tissue within the normal limits using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-XRCC4mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)





Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Adenocarcinoma of endometrium tissue using anti-XRCC4mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)







Immunohistochemical staining of paraffinembedded Carcinoma of prostate tissue using anti-XRCC4mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)

Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-XRCC4mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500650, Dilution 1:50)