

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA500773**

### DCXR Mouse Monoclonal Antibody [Clone ID: OTI7D11]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI7D11
Applications: IHC, WB

**Recommend Dilution:** WB 1:2000, IHC 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DCXR (NP\_057370) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.51 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 25.9 kDa

**Gene Name:** dicarbonyl and L-xylulose reductase

Database Link: NP 057370 Entrez Gene 51181 Human

Background: DCXR is an enzyme that has both diacetyl reductase (EC 1.1.1.5) and L-xylulose reductase (EC

1.1.1.10) activities

Synonyms: DCR; HCR2; HCRII; KIDCR; P34H; PNTSU; SDR20C1; XR


**Protein Families:** Druggable Genome

**Protein Pathways:** Metabolic pathways, Pentose and glucuronate interconversions






## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DCXR ([RC200023], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DCXR. Positive lysates [LY402535] (100ug) and [LC402535] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Carcinoma of liver tissue using anti-DCXRmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500773, Dilution 1:50)