

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500804

hnRNP F (HNRNPF) Mouse Monoclonal Antibody [Clone ID: OTI5F5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5F5

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:50, IF 1:50~100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HNRNPF (NP_004957) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 45.7 kDa

Gene Name: heterogeneous nuclear ribonucleoprotein F

Database Link: NP 004957 Entrez Gene 3185 Human

Background: This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear

ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNAs which have guanosine-rich sequences. This protein is very similar

to the family member hnRPH. Multiple alternatively spliced variants, encoding the same


protein, have been identified.

Synonyms: HNRPF; mcs94-1; OK; SW-cl.23







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HNRNPF ([RC205100], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HNRNPF. Positive lysates [LY401544] (100ug) and [LC401544] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 10 Human tissue by using anti-HNRNPF monoclonal antibody at 1:1000 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon;10: spleen).




Equivalent amounts of cell lysates (10 ug per lane) of wild-type 293T cells (WT, Cat# LC810293T) and HNRNPF-Knockout 293T cells (KO, Cat# [LC811367]) were separated by SDS-PAGE and immunoblotted with anti-HNRNPF monoclonal antibody TA500804, (1:500). Then the blotted membrane was stripped and reprobed with anti-HSP90AB1 antibody ([TA500494]) as a loading control.






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-HNRNPF mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500804)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-HNRNPF mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500804)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-HNRNPF mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500804)






Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-HNRNPF mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500804)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-HNRNPF mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500804)




Anti-HNRNPF mouse monoclonal antibody (TA500804) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HNRNPF ([RC205100]).






Immunofluorescent staining of HepG2 cells using anti-HNRNPF mouse monoclonal antibody (TA500804).



Flow cytometric Analysis of Hela cells, using anti-HNRNPF antibody (TA500804), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-HNRNPF antibody (TA500804), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).