

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500814

AKT2 Mouse Monoclonal Antibody [Clone ID: OTI4H7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4H7

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:1000, IHC 1:50, IF 1:50~100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human AKT2 (NP_001617) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.56 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 55.6 kDa

Gene Name: AKT serine/threonine kinase 2

Database Link: NP 001617 Entrez Gene 208 Human

Background: This gene is a putative oncogene encoding a protein belonging to a subfamily of

serine/threonine kinases containing SH2-like (Src homology 2-like) domains. The gene was shown to be amplified and overexpressed in 2 of 8 ovarian carcinoma cell lines and 2 of 15 primary ovarian tumors. Overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable

of phophorylating several known proteins. [provided by RefSeq]

Synonyms: HIHGHH; PKBB; PKBBETA; PRKBB; RAC-BETA

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase







Protein Pathways:

Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Tight junction, Toll-like receptor signaling pathway, VEGF signaling pathway

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AKT2 ([RC217733], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AKT2. Positive lysates [LY419836] (100ug) and [LC419836] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-AKT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500814)






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-AKT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500814)




Flow cytometric Analysis of Jurkat cells, using anti-AKT2 antibody (TA500814), (Red), compared to a nonspecific negative control antibody, (Blue).




Flow cytometric Analysis of Hela cells, using anti-AKT2 antibody (TA500814), (Red), compared to a nonspecific negative control antibody, (Blue).






Anti-AKT2 mouse monoclonal antibody (TA500814) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY AKT2 ([RC217733]).



Immunofluorescent staining of HepG2 cells using anti-AKT2 mouse monoclonal antibody (TA500814).



Western blot analysis of extracts (10ug) from 3 different cell lines by using anti-AKT2 monoclonal antibodyat 1:200 dilution.