

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500855

GRASP1 (GRIPAP1) Mouse Monoclonal Antibody [Clone ID: OTI2C11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2C11
Applications: FC, IF, WB

Recommend Dilution: WB 1:1000, IF 1:100, Flow 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human GRIPAP1 (NP_064522) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 96.0 kDa

Gene Name: GRIP1 associated protein 1

Database Link: NP 064522 Entrez Gene 56850 Human

Background: This gene encodes a guanine nucleotide exchange factor for the Ras family of small G

proteins (RasGEF). In brain studies, the encoded protein was found with the GRIP/AMPA receptor complex. Multiple alternatively spliced transcript variants have been described that encode different protein isoforms; however, the full-length nature and biological validity of all

of these variants have not been determined


Synonyms: GRASP-1

Protein Families: Druggable Genome







Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GRIPAP1 ([RC210936], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GRIPAP1. Positive lysates [LY402754] (100ug) and [LC402754] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (10ug) from 5 different cell lines by using anti-GRIPAP1 monoclonal antibody (1:200).



Anti-GRIPAP1 mouse monoclonal antibody (TA500855) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY GRIPAP1 ([RC210936]).



HEK293T cells transfected with either pCMV6-ENTRY GRIPAP1 ([RC210936]) (Red) or empty vector control plasmid (Blue) were immunostained with anti-GRIPAP1 mouse monoclonal (TA500855), and then analyzed by flow cytometry.