

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500880

XRCC1 Mouse Monoclonal Antibody [Clone ID: OTI2F8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F8

Applications: IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human XRCC1 (NP_006288) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.64 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 69.5 kDa

Gene Name: X-ray repair cross complementing 1

Database Link: NP 006288 Entrez Gene 7515 Human

Background: The protein encoded by this gene is involved in the efficient repair of DNA single-strand

breaks formed by exposure to ionizing radiation and alkylating agents. This protein interacts with DNA ligase III, polymerase beta and poly (ADP-ribose) polymerase to participate in the base excision repair pathway. It may play a role in DNA processing during meiogenesis and recombination in germ cells. A rare microsatellite polymorphism in this gene is associated

with cancer in patients of varying radiosensitivity.


Synonyms: RCC

Protein Families: Druggable Genome
Protein Pathways: Base excision repair






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY XRCC1 ([RC204952], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-XRCC1 (1:1000).




Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-XRCC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-XRCC1 mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-XRCC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)




Immunohistochemical staining of paraffinembedded lung tissue within the normal limits using anti-XRCC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-XRCC1 mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)






Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-XRCC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)





Immunohistochemical staining of paraffinembedded prostate tissue within the normal limits using anti-XRCC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of prostate tissue using anti-XRCC1 mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)







Immunohistochemical staining of paraffinembedded bladder tissue within the normal limits using anti-XRCC1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500880, Dilution 1:50)

Anti-XRCC1 mouse monoclonal antibody (TA500880) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY XRCC1 ([RC204952]).